

SuperBond Instant - Self Cure

- Self Cure bonds attachments to impacted teeth in 10 seconds.
- Kit includes 1.2cc Activator syringe and 1.2 cc Adhesive syringe.
- 10 applicator needles.
- Instructions included.

Cat. No. 557-764

1 kit - \$ 65.00 each
3 kits - \$ 60.00 each
5 kits - \$ 55.00 each



SuperBond Kit Cat. No. 557-754

COMPONENTS:

SuperBond Instant Activator only - Cat. No. 557-765 - \$35.00 each

SuperBond Instant Adhesive only - Cat.No. 557-766 - \$35.00 each

CuspidBond - Light Cure

- Light Cure bonds attachments to impacted teeth in 10 seconds.
- Kit includes 1.2cc Primer syringe and 1.2 cc Adhesive syringe.
- 10 applicator needles.
- Instructions included.

Cat. No. 557-754

1 kit - \$ 65.00 each
3 kits - \$ 60.00 each
5 kits - \$ 55.00 each



CuspidBond Kit Cat. No. 557-754

COMPONENTS:

CuspidBond Primer only - Cat. No. 557-755 - \$35.00 each

CuspidBond Adhesive only - Cat.No. 557-756 - \$35.00 each

Gold Chain and Gold Wire

- Gold Chain for cuspid traction.
- 12 " (30 cm) length.

Cat. No. 248-040

45.00 each



Gold Chain Cat. No. 248-040

- Gold Wire for tying attachments.
- 12 " (30 cm) length.

Cat. No. 248-045

10.00 each

Disposable Tips for Etch

- Perfectly fit onto 12 gm etch syringe to provide exact placement of etch.

Cat. No. 557-242

50 tips - \$ 8.00 each

100 tips - \$15.00 each



Autoclavable Lip Retractor


- Retracts lips and cheeks for bonding procedures.
- Soft flexible plastic for patient comfort.
- Autoclavable in paper bag.



Standard Expander
Cat. No. 822-444

Mini Expander
Cat. No. 822-666

1 to 4 expanders - \$ 6.00 each
5 to 9 expanders - \$ 5.40 each
10 + expanders - \$ 5.00 each



Coordinated use of Lip Retractor and Tongue Block

Tongue Block

- Prevents the tongue from interfering when etching and bonding.
- Bite block provides for patient's comfort during long procedures.
- Four assorted sizes per pack.

Cat. No. 822-111

1 to 4 packs - \$ 7.00 each
5 + packs - \$ 6.00 each



12 gm Green Syringe Etch

- For etching - apply to tooth enamel for 30 seconds.
- For hemostasis:
 1. Apply copiously to tooth and all hard and soft tissues for 30 seconds.
 2. Suction all etch away - do not wash.
 3. Reapply etch to entire operative field until edges of etch remain green (not black where bleeding continues).
 4. Suction, wash and meticulously dry.

CAT. NO. 557-618



1 to 9 syringes - \$ 6.00 each

10 + syringes - \$ 5.50 each

Mini Bonding Tweezers

- Easily picks up and securely holds all attachments.
- 5 different styles cover all placement situations.


1 to 4 - \$ 22.00 each

5 to 9 - \$ 20.00 each

10 + - \$ 18.00 each



Fine tips securely hold the tiniest attachment



Short - Cat. No. 878-020



Long - Cat. No. 878-021



Angled 30° - Cat. No. 878-022



Curved 45° - Cat. No. 878-023



Bayonet - Cat. No. 878-024

SPECIAL

Cat No 878-030

Five
Tweezers

One
of
each
style

\$ 90.00