

OralSurgerySource

Always the Best Products, Always the Best Price

13343 SHERMAN WAY
N. HOLLYWOOD, CA 91605 USA
TOLL FREE (888) 522 • 7223
INTERNATIONAL (818) 982 • 9445
24 HOUR FAX (818) 982 • 9501
e-mail OrthoSource@aol.com

CATALOG 47

ORTHODONTIC IMPLANTS - A MUST HAVE MODEL!

Demonstrates some possible ortho implant placements:

- Retraction of anterior teeth.
- Distalization without taxing anchorage.
- Anterior movement of posterior teeth
- Implant as molar anchorage to distalize bicuspid.
- Atraumatic placement of implants.

Cat. No. 840-125 \$ 255.00

**Spectacular Articulated
Orthodontic Implant Model
presents the advantages of
Implants in Orthodontics**



Implant Zygoma Base to Distalize Cuspid

Implant Zygoma Base to Distalize Cuspid




R
I
G
H
T

S
I
D
E



Ortho Implant Model - Anterior View



L
E
F
T

S
I
D
E

Implant Buccal Plate of Mandible to Mesialize Molar

Implant Posterior Mandible to Distalize Bicuspid

Gold Mesh disks - Perfect for Light and self cure Bonding!



Gold Mesh Disk Cat. No. 248-065

as low as

\$ 10⁰⁰
each

See page 3 for quantity pricing



CONTROL: Tips of Mini Tweezer fit into eyelet for complete control of attachment. Swivel allows chain to rotate full



ADHESIVE: With open Gold Mesh you can see resin locked into attachment guaranteeing bond to cuspid.



SHAPE: Wire cutter can trim Gold Mesh to fit any required surface. Mesh is easily bent to adapt to contour of any surface.



PLACEMENT: With tweezer control, ideal attachment shape, and adhesive in mesh, your bonding will be a success.

We will help you get into Light Cure CUSPID Bonding!

Light Cure Unit
LIST \$399

Light Cure Cuspid

- Rapid and Complete CurinG in 20 seconds.
- Autoclavable light guide turns 360°.
- Beep tone for curing times.
- Fan prevents bulb burn out.
- Available 110 V or 220 V



+



- Bonds attachment in Wet or Dry Field
- Rapid and Complete CurinG in 20 seconds through mesh.

**SUPER SPE-
CIAL
LIGHT & Resin
list \$464**